



SMALL EQUIPMENT

Autoclaves

MELAG

Vacuklav 31B "Class B" 18 Litres with printer **\$11,285.00**



Fractionated pre vacuum - removes air from the chamber
 Universal Program - Full Cycle 50 minutes. Innovative water system - single-use water principle.
 Post-Vacuum for best drying results. 4 Sterilization programs and 2 test programs. Easy to operate
 Dimensions - Outside D57cm x W42.5cm x H48.5cm. Inside 25cm x 35cm depth, 18 litres. Tray D29cm x W19cm x H2cm
 12 Months Warranty on Parts & Labour.

Vacuklav 23B "Class B" 22 Litres with printer (50 min) **\$11,650.00**

Vacuklav 41B 18 litres (34 min)

Chamber depth 35cm, includes 3 trays and on-board logger **\$13,800.00**

Vacuklav 43B 22 litres (34 min)

Chamber depth 45cm, includes 3 trays and on-board logger or to your networked PC **\$14,520.00**



 **Prestige Medical**



C3 Advance
 Faster Vacuum cycles
 6 tray chamber
 2 years unlimited usage warranty
 Data logger or PC link
 Manufactured in the UK

Superior quality Class B
 Large 16 or 22 litres
 Clean water system
 Optimal silent printer
 Compact design

12 Months Warranty on Parts & Labour.

16 litres with thermal printer **\$8,995.00**

16 litres with data logger **\$9,295.00**

22 litres with thermal printer **\$9,995.00**

22 litres with data logger **\$10,295.00**



MES

22 litres B & S Class with printer **\$8,723.00**

Chamber size 250mm x 450mm

2 year warranty

Three times pre vacuum processing vacuum dryness = .2% humidity.

Its universal B cycle sterilizes wrapped, unwrapped hollow or solid instruments.

Complete cycle in 49 minutes for wrapped instruments.

22 litres S Class with printer **\$7,147.00**

17 litres B & S Class with printer 36.5cm depth **\$5,999.90**

17 litres S Class with printer 36.5cm depth **\$5,490.00**



Ridley Dental Supplies are proud supporters of the Cancer Council and Oxfam

Order on line at www.ridleydentalsupplies.com.au