



NEW PRODUCTS

THERMIQUE Thermal Condenser Parkell, Inc

The Thermique thermal condenser only gets hot while its in the canal so users dont have to repeatedly heat their instrument in an open flame. The proportional feedback system maintains a precise temperature (+ 15C) at a precise depth within the canal for predictable gutta-percha flow and minimal generation of smoke and fumes. It also provides a large battery capacity and is compatible with vertical or lateral condensation.



\$1299.00

HANDPIECE LUBRICATION DEVICE

This machine can clean and oil the parallel shank or bending handpiece accurately and easily at any time.

1. Swing pinion can effectively clean and oil the dental handpiece from different angles.
2. Three dental handpieces can be oiled and cleaned at the same time.
3. Having completed the procedure of the cleaning and oiling, the remained oil can be cleared by pressing the air key.
4. Different handpieces can be oiled by choosing the short mode. The long mode, and the extra-long mode.

2x High Speed Connectors

1X Low Speed Connector

\$1199.95



CITIZEN Digital Ultrasonic Cleaner

- Features :
- Bigger tank.
 - Stronger transducer.
 - Digital timer&temp.
 - Heater Function.
 - Circuit Protector.
 - Wet-proof pcb.
 - Three separable wires.
 - Inside S.S.304 tank

Unit Size: 340x255x225mm

\$395.00



PROPHY MATE Polisher, Dental Descaler

Features:

1. Directly connecting to complex treatment unit, so it is very handy (2 hole & 4 hole connector only).
2. It is easy to replace the working head of Prophy-Mate. And it can be sterilized in condition of vacuum with temperature of 121.
3. Anti-resorption facility in san blasting gun can prevent sand blasting powder to be back-adsorbed into complex treatment unit.

\$395.00

